

LIFENET® System

Technical Specifications for the US

This information is being provided to help our customers configure their networks and prepare for installation of destination clients within their infrastructure.

CONNECTIVITY

LIFENET Applications make connections to the LIFENET System server by making outbound requests using TCP and TLS protocols using remote port 443. The destination applications make requests to the server using the SOAP data method to determine if new data exists for the destination application. The applications then pull the data to the destination and present it to the user in different methods depending on the type of destination and its configuration.

NETWORK CONFIGURATION

Internet access is required for the LIFENET System to be able to function properly. Network security systems may need to be configured to properly allow the LIFENET destination applications to communicate to the LIFENET System. There are many methods to secure a network that is connected to the internet. In order to help our customers successfully implement the LIFENET System in their networks, the following information should be considered.

For all solutions, it may be necessary to know details about where the LIFENET System is on the Internet. The FQDNs for the LIFENET System are <https://www.lifenetsystems.com>. There are two IP subnets for the LIFENET System. They are:

1. 209.67.69.176 255.255.255.240 or a 28 bit mask
2. 209.67.72.32 255.255.255.224 or a 27 bit mask

The system is co-located at a CenturyLink™ facility.

LIFENET SYSTEM PROPERTIES

Protocol	HTTPS
Connection	TLS 1.0, 1.1, 1.2; 128 bit or higher encryption; RSA 2048 bit exchange.
Address/URL	https://www.lifenetsystems.com/

WORKSTATION REQUIREMENTS

LIFENET System destination applications can be installed on existing workstations at your facility by users with administrative rights to that workstation. The workstations will need to have **.Net Framework 4.5** Microsoft Service Pack requirements will vary depending on the destination client being installed. Operating systems requirements will also vary depending on the destination client being installed. You may request this specific information from your LIFENET System sales or implementation representative. Clients support Microsoft® Windows® Vista Operating System or higher.

OTHER

Email filters may need to be set to allow emails originating from web@lifenetsystems.com.

For further information, please contact Physio-Control at 800.442.1142 (U.S.), 800.895.5896 (Canada) or visit our website at www.physio-control.com



Physio-Control Headquarters

11811 Willows Road NE
Redmond, WA 98052
www.physio-control.com

Customer Support

P. O. Box 97006
Redmond, WA 98073
Toll free 800 442 1142
Fax 800 426 8049