

**McGRATH**<sup>TM</sup>  
MAC **EMS**

Video Laryngoscope

**McGRATH™ MAC EMS**

Video Laryngoscope






**Thank you for choosing Physio-Control  
as your partner in helping save lives and  
improve patient care.**

The McGRATH™ MAC EMS video laryngoscope is shipped with one battery and the Instructions for Use. Disposable blades are shipped separately. In this catalog you will find an overview of blades and additional accessories available for the McGRATH™ MAC EMS. If you have any questions about accessories or disposables, please contact your Physio-Control representative.

## Disposable Laryngoscope Blades

---


McGRATH™ MAC Disposable Laryngoscope Blades are high performance, low-cost, single-use sterile blades. They are used for routine as well as many difficult airways, from pediatric to adult patients. They have a Macintosh-like curvature and are therefore easy to use with limited training. McGRATH™ MAC Blades are made of robust optical polymer and blade lenses are coated with hydrophilic optical surface treatment.



### McGRATH™ MAC 2 Disposable Laryngoscope Blade

Designed for pediatric patients older than eight weeks or at least 4.5 kg. Packed individually and shipped in boxes of 10.


**Box of 10 11996-000414**



### McGRATH™ MAC 3 Disposable Laryngoscope Blade

Designed for adult patients. Packed individually and shipped in boxes of 10.

**Box of 10 11996-000415**




### McGRATH™ MAC 4 Disposable Laryngoscope Blade

Designed for large adult patients. Packed individually and shipped in boxes of 10.

**Box of 10 11996-000416**

## McGRATH™ X blade™

The McGRATH™ X blade™ is a sterile disposable single-use blade that extends the range of proprietary disposable laryngoscope blades for the McGRATH™ MAC EMS system. The X blade™ contains a sweeping acute curvature for extreme airways, complementing the existing range of McGRATH™ MAC Disposable Laryngoscope Blades which are utilized in routine to difficult airways.



## Replacement Batteries



The McGRATH™ MAC EMS utilizes 3.6V Lithium power cells, which provide up to 250 minutes of usage. A representation of minutes remaining is displayed on screen. The on/off power switch is integrated into the battery pack, so when a new battery is fitted, the on/off function is renewed, keeping the product young.

### McGRATH™ 3.6V EMS Battery

The McGRATH™ MAC EMS is shipped with one 3.6V EMS Battery. Replacement batteries are shipped individually.  
**11996-000394**



# McGRATH™ MAC EMS

Video Laryngoscope





Physio-Control is now part of Stryker.

For further information please contact your local Physio-Control representative or visit our website at [www.physio-control.com](http://www.physio-control.com).

---

**Physio-Control Headquarters**

11811 Willows Road NE  
Redmond, WA 98052 USA  
Tel 425 867 4000  
[www.physio-control.com](http://www.physio-control.com)

**Physio-Control Operations  
Netherlands B.V.**

Galjoenweg 68  
6222 NV Maastricht  
The Netherlands  
Tel +31 (0)43 3620008  
Fax +31 (0)43 3632001

**Physio-Control  
UK Sales Ltd**

35 Great St Helen's  
London EC3A 6AP  
United Kingdom

**Physio-Control  
Australia Pty Ltd**

8 Herbert Street  
St Leonards NSW 2065  
Australia  
Toll Free Tel 1800 987 982  
Toll Free Fax 1800 890 892